





## **Description**

- High-quality treatment of all seed types using the rotor / stator principle
- Flexible treater, coating and palleting processes
- Perfect degree of treatment thanks to batch treatment technology
- Particularly suitable for seed breeding and application trials
- Perfect treatment results from 50 ml 2000 ml seed content
- Optimum grain-to-grain distributor
- Homogeneous grain coating
- Seed protection through air cushioning during treatment operation
- Fully visible outlet channel with excellent cleaning openings
- Spring-loaded discharge flap with on-hand operation
- Transportable and can be used anywhere (230 V and compressed air connection required)





## **Functions**

- Mixing cylinder with special mixing blades
- Infinitely variable speed adjustment of mixing cone
- Spray disk for fine spraying of the liquids (without nozzles)
- Treatment insertion by means of syringe, dispensette or peristaltic pump
- Very simple operation via manual switch
- Various discharge funnels, the right funnel for every bag
- Openings in the lid for aspiration and manual addition of treatment

## Optional:

- Storage container for seeds (approx. 1500 ml capacity according to customer requirement) with bottle flap
  - Labaratory peristaltic pump for standard treatment process
  - Completely closed mixing blades